

# Nelaton Catheters Standard

## DESCRIPTION

Our Nelaton Catheters are intended for intermittent bladder catheterisation. They are made from medical grade thermo-sensitive plastic and offer a perfectly finished distal end for smooth insertion. The blue coloured holding sleeve assists to maintain to maintain aseptic technique during insertion.

## FEATURES & BENEFITS

- Thermo-sensitive PVC
- Heat polished eyelets
- Smooth round tip
- Colour coded connectors
- Available for Paediatric, Female and Male patients
- Blue coloured holding sleeve to maintain aseptic technique during insertion


## PRODUCT DATA

- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Shelf life: 5 years
- ARTG number: 210992
- Materials: PVC
- Holding sleeve material: PE

## ISO STANDARDS

- ISO11607-1:2006 Packaging for terminally sterilised medical devices—Part 1: Requirements for materials, sterile barrier systems and packaging systems
- EN866:1999 Biology System for Testing Steriliser and Sterilisation Process
- ISO10993-10:2002 Biological evaluation of the medical devices-Part 10: Tests for irritation and sensitisation

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Tip
2	Side Hole
3	Tube
5	Blue Sleeve
4	Funnel - Colour Coded

# Nelaton Catheters Standard

Nelaton Catheters (Standard)					STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
<b>UR040000*</b>	Standard Nelaton Catheter Paediatric	30cm	6Fr	Light Green	50
<b>UR040001*</b>	Standard Nelaton Catheter Paediatric	30cm	8Fr	Light Blue	50
<b>UR040002</b>	Standard Nelaton Catheter Paediatric	30cm	10Fr	Black	50
<b>UR041000</b>	Standard Nelaton Catheter Female	20cm	8Fr	Light Blue	50
<b>UR041001</b>	Standard Nelaton Catheter Female	20cm	10Fr	Black	50
<b>UR041002</b>	Standard Nelaton Catheter Female	20cm	12Fr	White	50
<b>UR041003</b>	Standard Nelaton Catheter Female	20cm	14Fr	Green	50
<b>UR041004</b>	Standard Nelaton Catheter Female	20cm	16Fr	Orange	50
<b>UR041005</b>	Standard Nelaton Catheter Female	20cm	18Fr	Red	50
<b>UR041006</b>	Standard Nelaton Catheter Female	20cm	20Fr	Yellow	50
<b>UR042000</b>	Standard Nelaton Catheter Male	40cm	8Fr	Light BLue	50
<b>UR042001</b>	Standard Nelaton Catheter Male	40cm	10Fr	Black	50
<b>UR042002</b>	Standard Nelaton Catheter Male	40cm	12Fr	White	50
<b>UR042003</b>	Standard Nelaton Catheter Male	40cm	14Fr	Green	50
<b>UR042004</b>	Standard Nelaton Catheter Male	40cm	16Fr	Orange	50
<b>UR042005</b>	Standard Nelaton Catheter Male	40cm	18Fr	Red	50
<b>UR042006*</b>	Standard Nelaton Catheter Male	40cm	20Fr	Yellow	50

\* This product is currently being sold without the blue sleeve but will transition to blue sleeve upon depletion of existing stock.