



## Artemis Procedure Face Mask

The PrimeOn Artemis Procedure Face Mask provides excellent protection while remaining light and comfortable for everyday use. Resistant against fluid splashes, this latex-free design provides comfort and confidence for the wearer. Elastic earloops ensure the mask is easy to don and remove, yet will stay firmly in place whilst performing your daily tasks.



Middle Layer  
Meltblown filter

Inner Layer  
Fluid repellent to prevent bacterial growth

Outer Layer  
Additional fluid repellent material

### PrimeOn Artemis with Earloops

TGA Registration	ARTG 290961
Fixation	Earloops
Number of Fabric Layers	3
Anti-fog Foam	No
Latex Free	Yes
Colour	White
Fluid Resistance, mmHg	120
Bacterial Filtration Efficiency (BFE), %	≥98
Differential Pressure (ΔP), mm H <sub>2</sub> O/cm <sup>2</sup>	< 5.0

REORDER CODE  
FPM21UN

LEVEL BARRIER CLASSIFICATION  
Level 2

RECOMMENDED USE  
Where minimal blood or bodily fluid droplet exposure may occur

PACKAGING  
50 masks per box  
6 boxes per carton

STANDARDS  
AS4381-2015, ASTM F2100

MANUFACTURING ACCREDITATIONS  
ISO 9001, ISO 13485

A brand by



**Mun Australia**  
Toll free: 1800 456 837  
[munglobal.com.au](http://munglobal.com.au)



Protection Always On