



4-Compartment Procedure Tray



The GuardOn range of medical hollowware is made from sugarcane bagasse, a readily renewable resource and an environmentally friendly alternative to plastics because unlike plastic which take years to break down in landfill, it is compostable and biodegradable.

The GuardOn 4-Compartment Procedure Tray is an ideal solution for organising and holding a variety of medical essentials thanks to its versatile layout of one large compartment and three smaller separate compartments, which can accommodate a variety of equipment, fluids, medicine as well as instruments to suit different procedures according to their workflow. This tray's thoughtful design not only enhances medical practitioners' efficiency but also contributes to maintaining a hygienic and organised medical environment.

GuardOn 4-Compartment Procedure Tray

Dimensions

266 x 204 x 30mm

Material

Sugarcane Bagasse

Features

Rigid Base, leak & oil proof,
puncture resistant

Colour

Natural

Surface Finish

Smooth

Liquid Capacity

Holds liquids up to 90°C
for at least 24 hours

Liquid Compatibility

Not compatible with
alcohol based liquids

REORDER CODE

EHD14PT

STANDARDS

EN 13432

RECOMMENDED USE

Multipurpose and versatile, may be
used as an anaesthetic tray.

PACKAGING

250 units per carton

MANUFACTURING ACCREDITATIONS

ISO 14001 Environment Management

A brand by



Mun Australia

Toll free: 1800 456 837

munglobal.com.au



Protection Always On