

Kidney Dish



The GuardOn range of medical hollowware is made from sugarcane bagasse, a readily renewable resource and an environmentally friendly alternative to plastics because unlike plastic which take years to break down in landfill, it is compostable and biodegradable.

The GuardOn Kidney dish is particularly well-suited for applications involving medical instruments, wound dressing materials, or wipes, conveniently shaped to contain solutions or bodily fluids. An essential item for clinics, emergency departments, theatres and treatment rooms that is now more sustainable while enhancing the efficiency and convenience of medical procedures and hygiene routines.

GuardOn Kidney Dish

Dimensions	248 x 125 x 49mm
Material	Sugarcane Bagasse
Features	Rigid Base, leak & oil proof, puncture resistant
Colour	Natural
Surface Finish	Smooth
Liquid Capacity	Holds liquids up to 90°C for at least 24 hours
Liquid Compatibility	Not compatible with alcohol based liquids

REORDER CODE
EHD11KD

STANDARDS
EN 13432

RECOMMENDED USE
Ideally suited for medical instruments, wound dressing material or wipes

PACKAGING
400 units per carton

MANUFACTURING ACCREDITATIONS
ISO 14001 Environment Management

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au



Protection Always On