




Athena Surgical Face Mask

The PrimeOn Athena Surgical Face Mask provides protection at a premium level for surgical tasks. Our anti-fog foam technology prevents condensation from escaping out of the top of the mask, stopping eyewear from fogging, and ensuring maximum visibility at all times. The PrimeOn Athena surgical face mask is resistant against fluid splashes and conforms to the Australian Standard. The tie design ensures comfort for different head shapes and sizes.



PrimeOn Athena with Ties	
TGA Registration	ARTG 290961
Fixation	Ties
Number of Fabric Layers	4
Anti-fog Foam	Yes
Latex Free	Yes
Colour	Blue
Fluid Resistance, mmHg	160
Bacterial Filtration Efficiency (BFE), %	> 99
Differential Pressure (ΔP), mm H ₂ O/cm ²	< 5.0

REORDER CODE FPM30UN

LEVEL BARRIER CLASSIFICATION Level 3

RECOMMENDED USE Where the risk of exposure to blood or bodily fluid splash are high

PACKAGING 50 masks per box 6 boxes per carton

STANDARDS AS4381-2015, ASTM F2100

MANUFACTURING ACCREDITATIONS ISO 9001, ISO 13485





