



Item number | 1000021345



Pad, ABENA Let, Anatomic, blue colourcode, with SAP

- ✓ Solution for broad applications
- ✓ Support healthy skin balance
- ✓ Efficient core technology with Top-Dry Technology
- ✓ Discreet security with invisible leakage barriers
- ✓ Can be sterilized



Product Description

ABENA Let Anatomic is ABENA's anatomically shaped multipurpose pad with non-woven backsheet. The efficient core with Top-Dry Technology ensures skin integrity and optimal skin care by keeping the exposed area comfortable and dry even after multiple voids. The fully breathable pads are fitted with invisible lying leakage barriers for extra protection. The odour control system minimizes unpleasant smells. All products feature the Nordic Swan Ecolabel.











Specifications

Base name	Pad
Brand	ABENA
Sub-brand	Let
Size	Anatomic
Color	White
Features	Blue colourcode, with SAP
Material	ECF fluff, nonwoven, PP, resin, SAP
Length/depth	440 mm
Width	200 mm
Tolerance	+/- 10%
Absorption Rothwell	500 ml
Gender	Unisex
Certifications	UKCA, CE, MD, Nordic Swan Ecolabel
CE Class (Medical Devices)	Class I
Product or test standards	ISO 11948-1, ISO 13485, ISO 14001, ISO 9001
Directives, regulations and acts	MDR (EU) 2017/745
Safety Instructions and Warnings	Keep out of reach of children.
Because of the products composition it has an expiry date from the production date on:	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Dispose with household waste.
Packaging Disposal Instructions	Can be recycled or incinerated.



Instructions for use/application

Anatomically shaped multipurpose pad. Can be sterilized. It should always be used with special fixation underwear, like our ABENA Fix products. The fixation pants ensure a tight fit and reduce the risk of leakages.

Packaging data

Unit	Contains	Length	Width	Height	EAN
cll	270 pcs	73.6 cm	39.6 cm	31.2 cm	5713571001332
pck	30 pcs	22.6 cm	12 cm	29 cm	5713571001325
pcs	1 pcs				



Nordic Swan Ecolabel is an official Scandinavian environmental label. Products with the Nordic Swan Ecolabel meet stringent environmental requirements with regard to raw material, production, use, documentation, and safety.



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.



Odor control is a critical factor in ensuring discreteness for the user. ABENA's incontinence products are designed and engineered to ensure that minimal Odor arises from the product.



Unique Top-Dry Technology that provides a rapid absorption and ensures a dry surface even after multiple voids.



The product does not contain latex.



Drop system indicating absorption. The more drops, the higher the absorption. Ranges from 1 to 9 drops.





