

Intermittent Closed Catheterisation System

PRODUCT DATA SHEET



DESCRIPTION

Our Intermittent Closed Catheterisation System is sterile and pre-lubricated with Hydrophilic gel, so no additional preparations is needed. The closed system works like a traditional Nelaton Catheter and provides a simple catheterisation that reduces the risk of urinary tract infections.


FEATURES AND BENEFITS

- Sterile closed system
- Pre-lubricated
- Heat polished eyelets
- Colour coded
- Smooth round tip

PRODUCT DATA

- Bag size: 1500mL
- Side holes: 2
- Length: 42cm
- Material: PVC
- BPA: No
- Sterile EO
- DEHP free
- Latex free
- ARTG Number: 210992
- Shelf life: 3 Years

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	GAUGE	BAG CAPACITY	CARTON/QTY
UR050008	Intermittent Closed Catheterisation System	8Fr	1500mL	25pcs
UR050010	Intermittent Closed Catheterisation System	10Fr	1500mL	25pcs
UR050012	Intermittent Closed Catheterisation System	12Fr	1500mL	25pcs
UR050014	Intermittent Closed Catheterisation System	14Fr	1500mL	25pcs
UR050016	Intermittent Closed Catheterisation System	16Fr	1500mL	25pcs